

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Measurement Uncertainty: Summary Table for Immunoblot Assays

Measurement Uncertainty (MU) is calculated using internal quality control (IQC)

The raw data and calculations can be found at the following location:

<\\xggc-fsrv-04\GGC Biochemistry\South Glasgow\Immunology Common\COMMON\1IMM&NI\QUALITY\Quality Control\Uncertainty of Measurement>

Analyte	Neuronal antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)	
	Yo IQC	Titin IQC
Mean (x)	69.3	71.0
Number of measurements (n)	61	61
Estimated Standard Deviation (s)	12.75	8.24
Coefficient of Variance (%CV)	18.40%	11.61%
Coverage factor (k) to define a confidence level of 95%	2	2
Relative standard expanded uncertainty (U)	± 36.8%	± 23.2%
Uncertainty of measurement example (using IQC)	69.3 ± 25.5 signal intensity (43.8– 94.8 signal intensity) **Signal intensity range remains within positive qualitative group**	71.0 ± 16.5 signal intensity (54.5 – 87.5 signal intensity) **Signal intensity range remains within positive qualitative group**

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Analyte	Neuronal antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)	
	Hu IQC	CV2 IQC
Mean (x)	138.3	45.2
Number of measurements (n)	61	61
Estimated Standard Deviation (s)	3.31	6.95
Coefficient of Variance (%CV)	2.39%	15.38%
Coverage factor (k) to define a confidence level of 95%	2	2
Relative standard expanded uncertainty (U)	± 4.8%	± 30.8%
Uncertainty of measurement example (using IQC)	138.3 ± 6.6 signal intensity (131.7 – 145.0 signal intensity) **Signal intensity range remains within positive qualitative group**	45.2 ± 13.9 signal intensity (31.3 – 59.1 signal intensity) **Signal intensity range remains within positive qualitative group**

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Analyte	Neuronal antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)
	Amphiphysin IQC
Mean (x)	106.2
Number of measurements (n)	15
Estimated Standard Deviation (s)	18.18
Coefficient of Variance (%CV)	17.12%
Coverage factor (k) to define a confidence level of 95%	2
Relative standard expanded uncertainty (U)	± 34.23%
Uncertainty of measurement example (using IQC)	106.2 ± 36.4 signal intensity (69.8– 142.6 signal intensity) **Signal intensity range remains within positive qualitative group**

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Analyte	Liver antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)		
	M2 IQC	3E(BPO) IQC	Sp100 IQC
Mean (x)	94.5	128.9	60.8
Number of measurements (n)	18	18	18
Estimated Standard Deviation (s)	6.20	6.56	9.49
Coefficient of Variance (%CV)	6.56%	5.09%	15.61%
Coverage factor (k) to define a confidence level of 95%	2	2	2
Relative standard expanded uncertainty (U)	± 13.1%	± 10.2%	± 31.2%
Uncertainty of measurement example (using IQC)	94.5 ± 12.4 signal intensity (82.1 – 106.9 signal intensity)	128.9 ± 13.1 signal intensity (115.8 – 142.0 signal intensity)	60.8 ± 19.0 signal intensity (41.8 – 79.8 signal intensity)

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Analyte	Liver antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)		
	PML IQC	gp210 IQC	LKM IQC
Mean (x)	17.4	51.9	50.6
Number of measurements (n)	10	18	18
Estimated Standard Deviation (s)	3.31	13.16	8.09
Coefficient of Variance (%CV)	19.0%	25.36%	16.0%
Coverage factor (k) to define a confidence level of 95%	2	2	2
Relative standard expanded uncertainty (U)	± 38.0%	± 50.7%	± 32.0%
Uncertainty of measurement example (using IQC)	17.4 ± 6.6 signal intensity (10.8 – 24.0 signal intensity)	51.9 ± 26.3 signal intensity (25.6 – 78.2 signal intensity)	50.6 ± 16.2 signal intensity (34.4 – 66.7 signal intensity)

NHS Greater Glasgow & Clyde Immunology and Neuroimmunology		
QF_74	Measurement Uncertainty Summary Table for Immunoblot assays	Version: 3
Author: Lauren Hennessy	Authoriser: Sylvia Arthur / Fran Henriquez	Date of Issue: 27/06/25

Analyte	Liver antibodies using BlotMaster and EUROline scanner (assays reported qualitatively)		
	LC-1 IQC	SLA IQC	Ro52 IQC
Mean (x)	142.1	8.2	86.6
Number of measurements (n)	18	10	10
Estimated Standard Deviation (s)	3.48	2.20	4.48
Coefficient of Variance (%CV)	2.45%	26.84%	5.17%
Coverage factor (k) to define a confidence level of 95%	2	2	2
Relative standard expanded uncertainty (U)	± 4.9%	± 53.7%	± 10.3%
Uncertainty of measurement example (using IQC)	142 ± 7.0 signal intensity (135 – 149 signal intensity)	8.2 ± 4.4 signal intensity (3.8 – 12.6 signal intensity)	86.6 ± 9.0 signal intensity (77.6 – 95.6 signal intensity)

For the details of the calculation and MU protocol please refer to document QP_5: Measurement Uncertainty Protocol.