

<b>NHS Greater Glasgow &amp; Clyde Immunology and Neuroimmunology</b>		
<b>QF_4</b>	<b>Uncertainty of Measurement Summary Table for Optilite assays</b>	<b>Version: 4</b>
<b>Author: Lauren Hennessey</b>	<b>Authoriser: Moira Thomas</b>	<b>Date of Issue: 20/03/24</b>

### Uncertainty of Measurement: Summary table for Optilite assays

Uncertainty of measurement (UoM) is calculated using internal quality control (IQC)

The raw data and calculations can be found at the following location:

<\\xggc-fsrv-04\GGC Biochemistry\Immunology\1IMM&NI\Quality\Uncertainty of Measurement>

Analyte	C1 inhibitor concentration on Optilite	
	C1 Inhibitor Low IQC	C1 Inhibitor High IQC
Mean (x)	0.181	0.339
Number of measurements (n)	35	28
Estimated Standard Deviation (s)	0.008	0.014
Coefficient of Variance (%CV)	4.21%	4.26%
Coverage factor (k) to define a confidence level of 95%	2	2
<b>Relative standard expanded uncertainty (U)</b>	<b>± 8.4%</b>	<b>± 8.5%</b>
Uncertainty of measurement <b>example</b> (using IQC)	0.181 ± 0.015 g/L (0.166 – 0.196 g/L)	0.339 ± 0.029 g/L (0.310 – 0.368 g/L)

<b>NHS Greater Glasgow &amp; Clyde Immunology and Neuroimmunology</b>		
<b>QF_4</b>	<b>Uncertainty of Measurement Summary Table for Optilite assays</b>	<b>Version: 4</b>
<b>Author: Lauren Hennessey</b>	<b>Authoriser: Moira Thomas</b>	<b>Date of Issue: 20/03/24</b>

Analyte	<b>Serum free light chains (Kappa) on Optilite</b>	
	Kappa low positive IQC	Kappa high positive IQC
Mean (x)	44.0	384.7
Number of measurements (n)	228	161
Estimated Standard Deviation (s)	2.39	21.83
Coefficient of Variance (%CV)	5.43%	5.67%
Coverage factor (k) to define a confidence level of 95%	2	2
<b>Relative standard expanded uncertainty (U)</b>	<b>± 10.9%</b>	<b>± 11.3%</b>
Uncertainty of measurement <b>example</b> (using IQC)	44.0 ± 4.78 mg/L (39.2 – 48.8 mg/L)	384.7 ± 43.7 mg/L (341– 428 mg/L)

<b>NHS Greater Glasgow &amp; Clyde Immunology and Neuroimmunology</b>		
<b>QF_4</b>	<b>Uncertainty of Measurement Summary Table for Optilite assays</b>	<b>Version: 4</b>
<b>Author: Lauren Hennesy</b>	<b>Authoriser: Moira Thomas</b>	<b>Date of Issue: 20/03/24</b>

Analyte	<b>Serum free light chains (Lambda) on Optilite</b>	
	Lambda low positive IQC	Lambda high positive IQC
Mean (x)	43.4	434.4
Number of measurements (n)	100	100
Estimated Standard Deviation (s)	2.05	21.6
Coefficient of Variance (%CV)	4.73%	4.97%
Coverage factor (k) to define a confidence level of 95%	2	2
<b>Relative standard expanded uncertainty (U)</b>	<b>± 9.5%</b>	<b>± 9.9%</b>
Uncertainty of measurement <b>example</b> (using IQC)	43.4 ± 4.10 mg/L (39.3 – 47.5 mg/L)	434.4 ± 43.2 mg/L (391 – 478 mg/L)

For the details of the calculation and UoM protocol please refer to document QP\_5: Uncertainty of Measurement Protocol.