

Accessories available to you.

Additional accessories to assist you in every day situations such as home, at work, on the telephone or in busy environments.





Boost your hearing aid performance in challenging listening situations

With today's technology, hearing aids have made speech understanding in most situations enjoyable. But sometimes, hearing aids are simply not enough. Therefore Phonak has developed a wide range of wireless accessories to boost the performance of your hearing aids and keep you connected.

Whether listening in a noisy environment, at a distance or watching TV, Phonak offers hearing solutions for your everyday needs. Read more about our Roger™ portfolio and other hearing aid accessories.

I would like to improve...

	Changing the volume or program on my hearing aid	Hearing mainly one person speaking	Hearing in small meetings or with a couple of people talking around a table	Hearing the television
Roger On		✓	✓	✓
Roger Select		✓	✓	✓
Roger TableMic II				
Roger Clip-On Mic		✓		
Phonak PartnerMic		✓		
Phonak TV Connector				✓
Phonak RemoteControl	✓			

Listening to music clearly and without disturbing others	Hearing phone calls clearly from my mobile phone	Hearing a landline phone at home or at work	Hearing mainly one person in a busy pub or restaurant	Hearing in large meetings or several people talking around a table
✓		✓	✓	✓
✓		✓	✓	✓
				✓
			✓	
✓				

Roger On – Dedicated for conversations that occur in background noise and over distance



Experience the hearing performance with Roger On. With its ability to distinguish the direction of incoming speech and Pointing mode 2.0, users have the flexibility to focus on the conversations that truly matter.

It can be placed on a table, held in the hand or worn by a presenter. Personal settings are easily handled with the myRogerMic app.

Feature Overview:

- For group or one-to-one conversations
- Differentiate from which direction speech is coming
- Focus on specific conversations
- Operating time: 8 hours
- Operating range: 25 metres / 80 feet



Roger Select – Designed for noisy environments



The Roger Select is a state-of-the-art wireless microphone designed to help you clearly hear in noise. Placed in the middle of a table, it allows you to hear the conversation around you by automatically and seamlessly switching from one talker to another.

You can connect your TV, laptop and other multimedia devices via an audio input to hear directly.

Feature Overview:

- Fully automatic settings
- 360 degree hearing capabilities
- Audio input for connection to multimedia devices
- Manual touch keys to select directionality
- Operating time: 8 hours
- Operating range: 22 metres / 66 feet



TableMic II – Designed for complex meeting environments



This microphone is dedicated for working adults who participate in various meetings. It selects the person who's talking and switches automatically between the meeting participants.

Multiple Roger TableMic IIs can be paired up, making it ideal for large meetings. It can also transmit the sound of multimedia devices e.g. computers or TVs, to your hearing aids.

Feature Overview:

- MultiBeam technology for clearer directionality
- Remote control for muting and adjusting the range
- Audio input for connection to multimedia devices
- Operating time: 16 hours
- Operating range: 20 metres / 66 feet



Roger Neckloop – Universal receiver compatible with any hearing aid or cochlear implant with a T-coil




The Roger NeckLoop is worn around the neck by adults, students and children (above 8 years old) to give access to the complete Roger microphone transmitter range and the benefit of Roger technology in noisy situations and over distance.

Just turn on your Roger NeckLoop and use the display and battery level indicator to easily see the status of the receiver.

Feature Overview:

- Visual display for connection, battery level and volume control
- Rechargeable built in battery
- Detachable neckloop
- Operating time: 10 hours
- Operating range: 30 metres / 100 feet



Roger Clip-On Mic – A wireless partner microphone



This lightweight mini microphone is ideal for single conversations, plus it includes an audio input for listening to multimedia and TV connectivity.

It can also be used alongside other Roger Clip-On Mics in a microphone network.

Feature Overview:

- One-click connection – no programming needed
- Audio input for connection to multimedia devices
- Used as standalone or with other Roger microphones
- Operating time: 6 hours
- Operating range: 10 metres / 33 feet



TV Connector – Easy to use streaming device for TV and music



The TV Connector is a simple plug and play interface to TV's and other audio sources, supporting the direct connectivity of Phonak hearing aids. It's just like using wireless headphones

AirStream™ technology delivers top rated sound quality from any TV and stereo system¹. The TV Connector can simultaneously stream audio signals to an unlimited amount of connected hearing aids.

Feature Overview:

- Streaming from any TV and stereo system
- Easiest TV streaming device compared to other direct connectivity solutions¹
- Simultaneous streaming to an unlimited amount of hearing aids with one device
- Operating range: 15 metres / 50 feet



Phonak PartnerMic – Discreet and easy-to-use remote microphone




Phonak PartnerMic™ is a discreet and easy-to-use microphone optimal for one-to-one conversations. Attach it anywhere onto your partner's shirt, blouse or jacket without having to worry about the positioning.

Thanks to Phonak AirStream™ technology, speech is streamed to your hearing aids without an additional device.

Feature Overview:

- One-click connection – no programming needed
- Operating time: 6 hours
- Operating range: 10 metres / 33 feet



Phonak RemoteControl – Easy way to control your hearing aids



Pocket-sized and lanyard-equipped, this device was designed to provide maximum user convenience. It is the perfect accessory for adjusting the volume and programs on your hearing aids on-the-go.

Users can use this remote control to adjust the volume of their hearing aids using the up and down buttons. Furthermore, they can switch between programs by using the program button.

Feature Overview:

- Three large buttons for simple and convenient volume control and program change
- Compact and discreet
- Operating time: 12 months
- Operating range: 1 metre / 3.2 feet (typical)



¹ Legarth, S., Latzel, M. & Rodrigues, T. (2018). Media streaming: The sound quality wearers prefer, Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed October 17th, 2018.

You may be eligible to receive FREE accessories.

There are two fantastic schemes which offer support to individuals who experience hearing loss;

Access to Work - Individuals who are working, or looking for employment opportunities, may be able to receive funding support.

Disabled Students Allowance - Students entering Further Education may be eligible to receive funding to support them in their studies.

For more information, and to see if you're eligible to receive free support, visit www.hearinggrants.co.uk