

## Advice Request on Decontamination of Medical Devices

Complete and return to  
[ggc.infectioncontroldecontamination@ggc.scot.nhs.uk](mailto:ggc.infectioncontroldecontamination@ggc.scot.nhs.uk)

Date	Requesting Person	Contact Details
29/05/2017	Stewart Milne	24620
Name of device	Optiflow THRIVE	
Use (including type of patient and procedure)	Pre-oxygenation	
Reason for new device	New technique	
Manufacturer	Re-usable / Single Use?	
Fisher and Paykel	Waterbath replaced every 24 hours. Nasal cannula and filter are single patient use. In practice in GRI the whole kit will be single patient use	
Query	None. Been requested to do	
Details of manufacturers recommended decontamination	Dispose of single patient unit after every patient. Replace waterbath every 24 hours (in GRI the whole unit will be single use).	
Recommendation of Decontamination Group	The Decontamination Group class the whole unit as single-use therefore there are no issues with its use.	
Date	29/05/2017	