

Advice Request on Decontamination of Medical Devices

complete and return to
ggc.infectioncontroldcontamination@ggc.scot.nhs.uk

Date	Requesting Person	Contact Details
15/06/2017	John Gaffney	01475 504717
Name of device	Xario 200 Console and probe	
Use (including type of patient and procedure)	Ultrasound machine for TRUS	
Reason for new device	Service development	
Manufacturer		Re-usable / Single Use?
Toshiba		reusable
Query	Which cleaning product should be used to decontaminate ultrasound equipment?	
Details of manufacturers recommended decontamination	Cleanisept wipes	
Recommendation of Decontamination Group	The use of the Tristel 3-wipe process should be adopted for use with these probes. This has been carried out in other areas of NHSGGC. Contact details for the Tristel rep are LesleyMurray@tristel.com who can arrange teaching and provide order codes.	
Date	27/06/2017	

