



NHS GREATER GLASGOW & CLYDE


Continence Product Formulary

COMMUNITY ADULTS CORE




iD PRODUCT NAME	SIZE	WORKING ABSORBENCY	iD CODE	IMAGE	TENA EQUIVALENT
-----------------	------	--------------------	---------	-------	-----------------

iD Expert Light - Shaped Pad with adhesive strip for mild to moderate urinary incontinence
- must be worn with close fitting underwear

iD Expert Light Extra	35 x 12.5	150-200ml	5160040280		Comfort Mini Extra
iD Expert Light Maxi	40.5 x 15.5	250-350ml	5160050280		Comfort Mini Super

iD Expert Form - Shaped pad - for moderate to severe urinary incontinence
- Must be worn with the correct size stretch pants

iD Expert Form Normal	53x26/24	up to 500ml	5310155280		Comfort Normal
iD Expert Form Plus	64x33/31	500-625ml	5310260420UK		Comfort Plus
iD Expert Form Extra Plus	64x33/31	700-1000ml	5310270420UK		Comfort Super

iD Expert Comfort Fixation Pants - Measure waist to ensure correct size is prescribed
5 Pairs every 6 months. Please see manufacturers guide for fitting and washing instructions.

iD Expert Comfort Super - Small	40-80cm	N/A	5410100050		Tena Fix XS
iD Expert Comfort Super - Medium	60-100cm	N/A	5410200050		Tena Fix S
iD Expert Comfort Super - Large	80-120cm	N/A	5410300050		Tena Fix M
iD Expert Comfort Super - XLarge	100-160cm	N/A	5410400050		Tena Fix L
iD Expert Comfort Super - XXLLarge	140-180cm	N/A	5410500050		Tena Fix XL
iD Expert Comfort Super - XXXLarge	180-195cm	N/A	5410600030		Tena Fix XXL

If a patient requires a product outwith this formulary please email
NHSGGC.SphereBBSservice@ggc.scot.nhs.uk